


Estimation of Hardness of Water



Advantages of EDTA method

This method is preferable because of


- 1) Greater accuracy
- 2) Convenience
- 3) More rapid procedure




EDTA method



- Ethylene diamine tetra acetic acid disodium salt (EDTA disodium salt) is used as a strong complexing agent with Ca^{2+} and Mg^{2+} in hard water.
- The structure of EDTA disodium salt is:



- Initially, Ca^{2+} and Mg^{2+} are treated with Eriochrome black T (EBT) indicator using ammonia buffer (to maintain pH between 9-10) to get an unstable complex of Ca^{2+} and Mg^{2+} formed with EBT.






Estimation Water Hardness



On addition of EDTA, EBT gets replaced by EDTA since EDTA forms a stronger complex with the metal ions.

This is indicated by the formation of a *steel blue* coloured complex.



Experimental Procedure



1) Standardization of EDTA

- First EDTA Solution is standardized using standard hard water (1 mg/mL of CaCO_3 equivalents is prepared as standard hard water).
- For this, first known aliquot of Standard hard water is taken and 10-15 mL of ammonia buffer is added to bring the pH between 9-10.
- Then a few drops of EBT solution is added to form the unstable complex giving **wine red** colour.
- This solution is titrated with the EDTA solution till the solution turns to **steel blue** indicating the formation of stable EDTA-Metal ion complex.
- This volume of EDTA is noted as V_1 .

Experimental Procedure



2) Estimation of Total Hardness

- The above procedure is repeated with sample hard water of unknown hardness.
- Volume of EDTA is noted as V_2 .

3) Estimation of Permanent Hardness

- Then sample hard water of 250 mL is taken and evaporated to a volume of 50mL when the temporary hard salts settle down.
- The solution is filtered and washed thoroughly and made up again to 250mL.
- From this solution, 50 mL is pipetted out and titrated in similar manner as done with standard hard water.
- Volume of EDTA is noted as V_3 .

Experiment

1. Preparation of EDTA solution (0.01 M)

i.e. dissolve 3.72 g of sodium salt of EDTA crystals in 1 litre of distilled water

2. Preparation of standard Hard Water

a. dissolve 1.0 g of pure, dry CaCO_3 in small quantity of conc. HCl and evaporate the solution to dryness

b. Dissolve the residue in distilled water (1 litre)

1ml of this solution contains 1 mg of CaCO_3 equ. hardness

3. Preparation of buffer solution

Add 67.5 g of NH_4Cl to 570 ml of conc. Ammonia soln. and dilute with distilled water to 1 litre

4. Preparation of indicator

Dissolve 0.5 g of EB-T in 100 ml alcohol



5. Standardization of EDTA solution

Burette - EDTA solution (0.01 M)



Conical Flask - 50 ml of std. hard water + 10-15 ml of ammonia buffer + few drops of EB-T indicator

Volume of EDTA consumed = V_1 ml

a) Standardization of EDTA:

V_1 mL of EDTA is consumed by 50 mL of std. hard water

V_1 mL of EDTA = 50 mg of CaCO_3

$$1 \text{ mL of EDTA} = \left(\frac{50}{V_1} \right) \text{ mg of } \text{CaCO}_3$$



6. Estimation of total Hardness

Conical Flask - **50 ml of sample hard water** + **10-15 ml of ammonia buffer** + **few drops of EB-T indicator**

b) Total hardness:

Volume of EDTA consumed = V_2 ml

$$1 \text{ mL of EDTA} = \left(\frac{50}{V_1} \right) \text{ mg of CaCO}_3$$

$$\text{So, } V_2 \text{ mL of EDTA} = \left(\frac{50}{V_1} \right) \times V_2 \text{ mg of CaCO}_3$$

$$\text{Therefore, 1000 mL of sample hard water} = \left(\frac{50}{V_1} \right) \times \left(\frac{V_2}{50} \right) \times 1000 \text{ mg/L}$$

$$\text{i.e. Total hardness of sample hard water} = \left(\frac{V_2}{V_1} \right) \times 1000 \text{ mg of CaCO}_3 \text{ (ppm.)}$$



7. Estimation of permanent hardness

Sample water is boiled and filtered

Conical Flask – **50 mL of boiled & filtered sample hard water** **+** **10-15 mL of ammonia buffer** **+** **Few drops of EB-T Indicator**

c) Permanent hardness:

Volume of EDTA consumed = V_3 mL

$$1 \text{ mL of EDTA} = \left(\frac{50}{V_1} \right) \text{ mg of CaCO}_3$$

$$\text{So, } V_3 \text{ mL of EDTA} = \left(\frac{50}{V_1} \right) \times V_3 \text{ mg of CaCO}_3$$

$$\text{Therefore, 1000mL of sample hard water} = \left(\frac{50}{V_1} \right) \times \left(\frac{V_3}{50} \right) \times 1000 \text{ mg/L}$$

$$\text{i.e. Permanent hardness of sample hard water} = \left(\frac{V_3}{V_1} \right) \times 1000 \text{ mg of CaCO}_3 \text{ (ppm.)}$$

Temporary hardness calculation



Temporary hardness:

Temporary hardness = Total hardness – permanent hardness

$$= \left(\frac{V_2}{V_1} \times 1000 \right) - \left(\frac{V_3}{V_1} \times 1000 \right)$$

$$= \left(\frac{V_2 - V_3}{V_1} \times 1000 \right)$$