

Scanning Tunneling Microscope (STM)

A diagram consisting of three circles arranged in a triangle. The top-left circle is blue and contains the text 'QUANTUM TUNNELING EFFECT'. The top-right circle is red and contains the text 'Tunnelling Microscope'. The bottom circle is purple and contains the text 'STM'. The circles are connected by lines, forming a triangle.

**QUANTUM
TUNNELING
EFFECT**

**Tunnelling
Microscope**

STM


Why STM ?



- The electronic microscopes gives 'volume images' (penetration depth)
- In STM-no use of external particles
- Principle-Electrons tunneling between an atomically sharp tip and a surface

The STM combines three main concepts:


- ***Scanning***
- ***Tunneling***
- ***Tip-point probing***

STM System Components






- **In STM - Bias voltage is applied between a sharp conductive tip and a conductive sample.**
- **Sample is approached to a few angstroms from the tip, tunneling current occurs, that indicates proximity of the tip to the sample with very high accuracy.**
- **Scanning tunneling microscopy can be used to study conductive surfaces or thin nonconductive films and small objects deposited on conductive substrates.**




Operating modes

- Constant height mode
 - Current changes exponentially
 - Requires a smooth surface
- Constant current mode
 - Beware of insulators


Constant height mode



Constant Current mode



(a) Constant Current Mode



(b) CONSTANT CURRENT MODE




Fig. 1.60. Schematic illustration of the constant current mode of STM operation (Binnig and Rohrer, 1987; Hansma and Tersoff, 1987).

Advantages

- No damage to the sample
- Vertical resolution superior to SEM
- Spectroscopy of individual atoms
- Relatively Low Cost

Disadvantages

- Samples limited to conductors and semiconductors
- Limited Biological Applications: AFM
- Generally a difficult technique to perform


Figures of Merit

- Maximum Vertical Resolution: $.1 \text{ \AA}$
- Maximum Lateral Resolution: 1 \AA
- Maximum Field of View: $100 \text{ }\mu\text{m}$

Applications of STM

- Surface Structure: Compare to bulk structure
- Stuff Physicists Do: Semiconductor surface structure, Nanotechnology, Superconductors, etc.
- Metal-catalyzed reactions
- Spectroscopy of single atoms
- Limited biological applications: Atomic Force Microscopy
- Future Developments: Improve understanding of how electronic structure affects tunneling current, continue to develop STM offshoots

Scanning Tunneling Microscope (STM) 1981 Binnig and Rohrer




How an STM works ...



© Michael Schmid
Institut f. Allgemeine Physik
TU Wien 1997-2002

The STM can obtain images of conductive surfaces at an atomic scale of 0.2 nm, and also can be used to manipulate individual atoms, trigger chemical reactions, or reversibly produce ions by removing or adding individual electrons from atoms or molecules.

STM allows manipulation of individual atoms (1989)



Xenon atoms spell IBM on a nickel plate (IBM)



Iron atoms spell "Atom" on copper in Kanji characters. (IBM)