


# Dissolved Oxygen Determination



- Sample water with DO + Chemical Compound which takes oxygen from water and oxidizes


Concentration of oxidized  
chemical compound  $\propto$  DO

# DO Determination - Winkler's method



i) Oxygen present in sample water oxidizes  $\text{Mn}^{2+}$  to  $\text{Mn}^{4+}$

# DO Determination - Winkler's method



i)  $\text{Mn}^{4+}$  in sulphate form reacts with KI and liberates  $\text{I}_2$   
(Iodometric titration)

ii)  $[\text{DO}] \propto [\text{Mn}^{4+}] \propto [\text{I}_2]$


*(Sodium thiosulphate can be used to estimate the concentration of  $\text{I}_2$  with Starch Indicator)*

# Estimation of DO by Winkler's method



## Theory of Winkler's Method:


- Oxygen in the water sample oxidizes iodide ion ( $I^-$ ) to iodine ( $I_2$ ) quantitatively.
- The amount of iodine generated is then determined by titration with a standard thiosulfate ( $S_2O_3^{2-}$ ) solution.
- The endpoint is determined by using starch as a visual indicator.







# Estimation of DO by Winkler's method



## Titration of Iodine Solution



# Estimation of DO by Winkler's method



- The amount of oxygen can then be computed from the titre values