

Electron Diffraction -

Apparatus Required

1. Electron diffraction tube with stand.
2. High voltage power supply (upto 10 kV).
3. Connecting wires.
4. Plastic measuring scale.

OBJECTIVE:-

To calculate the interplanar spacing in graphite from the diffraction pattern.

BASIC INFORMATION:

In this experiment, electrons get transmitted through a very thin polycrystalline graphite sheet. The schematic sketch is shown in fig 1. Graphite has two independent lattice spacings (d_1 and d_2) and these are shown in Fig. 2. The two diffraction rings that will be seen at each voltage are due to these planes.

Applying the diffraction formula for first order, we get ,

$$\lambda = d \sin \theta.$$

where λ is the de-Broglie wavelength of e^- , d is interplanar spacing and θ is the angle of diffraction. Electrons are accelerated through a p.d. of V Volts and hence their de-Broglie wavelength is ,

$$\lambda = \frac{12.3}{\sqrt{V}} \text{ \AA}$$

From the geometry of Fig. 1 we have,

$$\sin \theta = \frac{R}{\sqrt{R^2 + L^2}}$$

Upon simplifying and using the fixed value of $L = 13.5 \text{ cm}$ and R is expressed in cm,

$$\sin \theta = \frac{1}{\left(1 + \left(\frac{13.5}{R}\right)^2\right)^{0.5}}$$

Interplanar spacing can be calculated from equation (1) by substituting equations (2) and (4) into it.

SAFETY GUIDELINES AND PRECAUTIONS:-

1. Never accelerate beyond 5KV.
2. Never touch any controls on the power supply other than the 'on-off' switch and voltage varying knob.
3. Never use force to measure the ring diameters. Keep a plastic scale very gently over the tube to measure the diameters. Metal scales are not allowed.
4. You are working with a very high energy source (>5KV) and hence touching any part of the entire setup other than what is mentioned in point 3 (just for the purpose of measurement) is prohibited. This is for the purpose of your own safety and the safety of the Lab.

Procedure:-

1. Set the accelerating voltage at 4 kV.
2. For the inner ring, measure the diameter ($2R_1$).
3. Fill up the radius (R_1) in the tabular column.
4. For the outer ring, measure the diameter ($2R_2$).
5. Fill up the radius (R_2) in the tabular column.
6. Calculate λ , $\sin\theta$ and d from the equations (2), (4) and (1) respectively and fill up the corresponding cells in the tabular column.
7. Repeat steps 2 to 6 for accelerating voltages 4.5 and 5 kV.
8. Calculate the average 'd' for both inner and outer rings.

~~Result~~




Fig:- Schematic diagram of experimental setup.




Fig:- Left (Two independent types of lattice planes in polycrystalline .

Tabular Column:-

Ring	V (kv)	$2R_1$ or $2R_2$ (cm)	R_1 or R_2 (cm)	λ (Å)	$\sin\theta$	d (Å)
Inner	4.0	2.4	1.2	0.1945	0.0885	2.138
	4.5	2.2	1.1	0.1834	0.0812	2.259
	5.0	2.0	1.0	0.1739	0.0739	2.353
Outer	4.0	4.2	2.1	0.1945	0.1537	1.265
	4.5	3.9	1.95	0.1834	0.1430	1.282
	5.0	3.7	1.85	0.1739	0.1358	1.280

SAMPLE CALCULATION:-

For Inner radii, $V = 4.0$ kV

$$2R_1 = 2.4 \text{ cm}$$

$$R_1 = \frac{2.4}{2} = 1.2 \text{ cm}$$

$$\lambda = \frac{12.3}{\sqrt{V}} \text{ Å} = \frac{12.3}{\sqrt{4000}} \text{ Å} = 0.1945 \text{ Å}$$

$$\sin\theta = \frac{1}{\left(1 + \left(\frac{13.5}{R}\right)^2\right)^{0.5}} = \frac{1}{\left(1 + \left(\frac{13.5}{1.2}\right)^2\right)^{0.5}} = 0.0885$$

$$d = \frac{\lambda}{\sin\theta} = \frac{0.1945}{0.0885} \text{ Å} = 2.198 \text{ Å}$$

For outer ring,

$$\text{Averaged } d = \frac{2.198 + 2.259 + 2.353}{3} = 2.27 \text{ Å}$$

For inner ring,

$$\text{Averaged } d = \frac{1.265 + 1.282 + 1.280}{3} = 1.276 \text{ Å}$$

Similarly, inner diameter and outer ones were calculated and

observed by mean that the interplanar spacing are $d_1 = 2.27 \text{ Å}$

$$d_2 = 1.276 \text{ Å}$$

Date: 2021/04/01.

Reg:- 20BDE0405

Date

Page

RESULT:-

The Interplanar spacings in graphite were measured as $d_1 = 0.227 \text{ nm}$ and $d_2 = 0.1276 \text{ nm}$.