

INTEGRATED OPTICS.

APPARATUS AVAILABLE:-

- Spectrometer
- Spirit Level
- Magnifying Glass
- Glass prism
- Sodium Vapour Lamp.

SLO:

✓ To determine the apex angle of given prism using a spectrometer.

Date:- 2021/03/04

PHY 1701 (Engineering Physics)

Lab Manual and Record
Reg No:- 20BDS0405




Fig:- Schematic diagram of a Spectrometer.




Fig:- Ray Diagram for Angle of Prism.

Date: - 2021/03/04

Tabulation: -

$$\text{Least Count} = 1' = \left(\frac{1}{60}\right)^\circ$$

Reading of reflected ray	Vernier A			Vernier B		
	MSR	VSR	Total	MSR	VSR	Total
Reflection from side AB (a)	293°	17'	293.28°	113°	13'	113.22
Reflection from side AC (b)	55°	12'	55.2° + (360°)	235°	25'	235.42
Difference between 'a' and 'b'	$\theta_1 =$		121.92°	$\theta_2 =$		122.2

~~Mean =~~

SAMPLE CALCULATIONS: -

For Vernier A, for reflection of ray from AB, MSR = 293°, VSR = 17'

$$\text{Total} = 293 + \frac{17}{60} = 293.28^\circ$$

for reflection of ray from AC, MSR = 55°, VSR = 12'

$$\text{Total} = 55 + \frac{12}{60} = 55.2^\circ + (360^\circ)$$

$$\theta_1 = (360 + 55.2) - 293.28 = 121.92^\circ$$

For Vernier B, for reflection of ray from AB, MSR = 113°, VSR = 13'

$$\text{Total} = 113 + \frac{13}{60} = 113.22^\circ$$

for reflection of ray from AC, MSR = 235°, VSR = 25'

$$\text{Total} = 235 + \frac{25}{60} = 235.42^\circ$$

$$\theta_2 = (235.42 - 113.22) = 122.2^\circ$$

$$\text{Mean}(\theta) = \frac{\theta_1 + \theta_2}{2} = \frac{122.2 + 121.92}{2} = 122.06^\circ$$

$$\text{Angle of prism (A)} = \frac{\theta}{2} = \frac{122.06}{2} = 61.03^\circ$$

PHY 1701 (Engineering Physics)

Lab Manual and Records

Date:- 2021/03/04

Reg No:- 20BDS0405

RESULT:-

The apex angle of the given equilateral prism is 61.05°