


Date:Reg. No.MONOCHROMATORS IN SOPHISTICATED INSTRUMENT**Apparatus Available:**

- Laser Source
- Grating
- Scale with measurements

**SLO:**

- ✓ To determine the wavelength of the given laser source using transmission diffraction grating method



$D$  = the distance from the grating to the screen.

$d$  = the spacing between every two lines (same thing as every two sources)

If there are ( $N$ ) lines per mm of the grating, then ( $d$ ), the space between every two adjacent lines or (every two adjacent sources) is

$$d = \frac{1}{N}$$

The diffraction grating formula for the principal maxima is:  $d \sin \theta = n \lambda$

where  $n$  is the order of diffraction ( $= 1, 2, 3, \dots$ ) and  $\theta$  angle of diffraction.

$$\lambda = \frac{d \sin \theta}{n} \text{ (meter)}$$

**Observation:**

Number of lines per meter on the grating is \_\_\_\_\_


Date:Reg. No.

Fig 2. Experimental Setup

**Tabular Column:**

n	S (cm)	2L (cm)	L (cm)	$\tan\theta=L/S$	$\theta=\tan^{-1}(L/S)$	$\sin\theta$	Mean	$\lambda$ (nm)
1	25							
	30							
	35							
2	25							
	30							
	35							
3	25							
	30							
	35							

**Date:****Reg. No.****Calculation:****Result:**

The wavelength of the laser source is found to be \_\_\_\_\_ (nm)

=====