

Date:Reg. No.QUALITY CHECK FOR SOFT DRINKS**Apparatus Required:**

- Travelling microscope
- Transparent liquid (water)
- Reading lens
- Glass beaker
- Pin
- Saw dust

SLO:

To determine the refractive index of the given transparent liquid using travelling microscope

To determine the refractive index of impure liquid

Formula:

The refractive index of liquid


$$\mu = \frac{\text{Real depth of the liquid}}{\text{Apparent depth of the liquid}} = \frac{(C - A)}{(C - B)} \quad (\text{No Units})$$

where

A is the microscopic reading when tip of the pin is focused directly

B is the microscopic reading when tip of the pin is focused through the liquid

C is the microscopic reading when saw dust sprinkled on the surface of the liquid is focused



Date:

Reg. No.

Least count of travelling microscope = _____													
Table 1	Volume of water in the beaker	Clear image of tip of the pin (Reading A)			Clear image of tip of the pin seen through the liquid (Reading B)			Clear image of the saw dust scattered on the surface of liquid. (Reading C)			C-A	C-B	μ
		MSR ()	VSR ()	OR ()	MSR ()	VSR ()	OR ()	MSR ()	VSR ()	OR ()			
	40 ml												
	60 ml												
Mean													
VSR = VSC \times LC ; Observed Reading = MSR + VSR													

Date:

Reg. No.

Result:

Refractive index of the given liquid (water) is found to be _____ (Units)
