

## HELIUM-NEON LASER

#### Dr Rajeshkumar Mohanraman

Assistant Professor Grade 1 School of Advanced Sciences VIT Vellore

## Main Components

Three main components of ANY lasers are

- (i) The active medium
- (ii) The pumping source
- (iii) The optical resonator
- The active medium acts as an amplifier for light waves
- For amplification, medium should be in a state of Population inversion
- Population inversion metastable levels lifetime is bit longer as compared to excited state
- The active medium placed inside an optical resonator acts as an oscillator
- A pair of mirrors + active medium optical resonators

## Helium – Neon laser

- ❖ First continuous laser developed by Ali Javan, W. Bennutt and D. Herriot in 1961
- Operation wavelength is 632.8 nm (red portion of visible spectrum).
- ❖ 4-level laser scheme.
- ❖ More directional and monochromatic than solid state lasers.
- Output is moderate compared with solid state lasers.
- ❖ Active medium is a mixture of He and Ne gas in 10:1 ratio (atomic percentage).
- ❖ Ne atoms are active centres for lasing action. He only helps in efficient excitation of Ne atoms.






Completely protected commercial He-Ne laser


He-Ne laser – without protective jacket

### **Construction of He-Ne laser**



- ✓ Set up consists of a discharge tube of length 80 cm and bore diameter of 1.5 cm.
- ✓ Gain medium of the laser is a mixture of He and Ne as the name suggests in ~10:1 ratio. It is contained at low pressure (an average 50 Pa/cm of cavity length) in a glass envelope.
- ✓ The pumping is provided electrically by creating an discharge. The electrical discharge is created by applying
  ~ 1 KV through an anode and cathode present at each end of the glass tube. The typical current value ranges from 5 100 mA for continuous mode operation.

## He-Ne Energy level diagram




 $He + e_1 \rightarrow He^* + e_2$ 

 $He^* + Ne \rightarrow He + Ne^*$ 

- ✓ When voltage is applied to the electrodes it **ionizes the gas, the electrons and ions** thus produced are accelerated towards anode and cathode respectively.
- ✓ Electrons acquire higher velocity due to its smaller mass when compared to the others. They transfer **K.E** to **He atoms through inelastic collision**.
- ✓ He atoms are readily excited by electrons impact because of its fairly light mass.
- ✓ Thus He atoms are excited to  $F_2$  and  $F_3$  states which lie at 19.81 and 20.61 eV respectively.
- ✓ These are **metastable states** and these atoms cannot return to ground state readily by spontaneous emission.

- ✓ These atoms return to ground state by transferring energy to Ne atom in the state which has identical energy. Such transfer is called resonant transfer of energy. (The direct excitation of Ne atoms is in-efficient compared to He)
- ✓ Neon energy levels  $E_6$  and  $E_4$  nearly coincide with  $F_3$  and  $F_2$  of Helium, so resonant transfer can occur.
- ✓ The additional energy 0.05 eV is provided by the K.E of the He atom.
- ✓ This energy exchange process occurs with high probability only because of the
  accidental near equality of the two excitation energies of the two levels in these
  atoms. Thus, the purpose of population inversion is fulfilled.
- $\clubsuit$  Ne atoms in the  $E_6$  level and  $E_4$  level emit a photon parallel to the axis of the tube.
- ❖ This photon travels through the gas mixture parallel to the axis of tube, it is reflected back and forth by the mirror ends until it stimulates an excited Ne atom and causes it to emit a photon with the stimulating photon.
- ❖ In reality neon energy levels E<sub>6</sub>, E<sub>5</sub>, E<sub>4</sub>, E<sub>3</sub>, E<sub>2</sub> are not single but a group of lines. Consequently several laser transitions are possible.
- ❖ Three main laser transitions are
  - 1.  $E_6$  to  $E_3$  generates laser beam of red colour at 6328Å
  - 2.  $E_4$  to  $E_3$  IR beam at wavelength of 1.15  $\mu m$
  - 3.  $E_6$  to  $E_5$  light in Far IR region at 3.39  $\mu$ m

#### He-Ne Laser Spectrum



# **Applications of He-Ne laser**

- The Narrow red beam of He-Ne laser is used in supermarkets to read bar codes.
- The He- Ne Laser is used in Holography in producing the 3D images of objects.
- He-Ne lasers have many industrial and scientific uses, and are often used in laboratory demonstrations of optics.

# He-Ne lasers uses transitions among the various excited electronic states of an atom.

Advantages	Disadvantages
Emits laser light in the visible portion of the spectrum	Low efficiency
High stability	Low gain
Low cost	Output power is small
Operates without damage at higher temperatures	

#### **CAUTION**



Helium-neon lasers are common in the introductory physics laboratories, but they can still be dangerous! According to Garmire, an unfocused 1-mW HeNe laser has a brightness equal to sunshine on a clear day (0.1 watt/cm²) and is just as dangerous to stare at directly.